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IT – CUFFIE PER EEG DISPOCAP MONOUSO

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Instruction for use

EN – DISPOSABLE DISPOCAP EEG HEADCAP

DISPOCAP differ from traditional headcap for its versatility and the advanced solutions that offer the user. The main advantages are:

- The net of the headcap is built in an elastic fabric mesh very large to provide optimum comfort to the patient. The holding of the net on the patient's head is made by a silicone rubber strap acting on the entire circumference of the head, so it is not necessary to use under-chin and thoracic belts.
- The electrodes may be of Tin, Silver, Silver/Silver Chloride or Silver/Silver Chloride sintered Ag/AgCl.
- The electrode and the cable are moulded with plastic material, being one piece only it does not allow the migration of gel towards the cable.
- The electrodes are divided into two screwed parts and may be easily mounted by the operator in the preferred position.
- The electrode looks like a sole body, with a central hole to inject the conductive paste. The net of the cap is not pressed on the holding of the electrode holder, that is why it does not come into direct contact with the gel; the net is put into the electrode and blocked with a ring;

The cap is supplied with two cables made of twisted wires of different colour with 10 electrodes and a single cable with electrode to be used as a reference.

All the individual wires have T.P. DIN 42802 Connectors.

Cables with different numbers of electrodes can be supplied upon request.

The cap is supplied in different sizes.

In order to obtain an optimal recording please take note of the following instructions:

1. Choose the cap most suitable to the patient's head.
2. Make sure that the cap fits the patient's head uniformly.
3. Fill with gel (**NEURGEL**) the syringe equipped with a blunted needle, insert the needle in the electrode hole and rub slightly the skin with a round movement.
4. Inject a sufficient quantity of gel to maintain a constant contact between the patient's skin and the electrode.
5. Make sure that the skin impedance is not too much high.
6. In case of high impedance on some electrodes repeat what is described at point 3.
7. A good impedance (lowest level possible) on the ground and reference electrodes offers optimal recordings conditions, for this reason it is advisable, if used a ground or reference ear electrode, to use electrodes made by the same material as the recordings electrodes.

In case of artefact: the cap could be too flabby, some residual dry gel could be on the electrode, the ground electrode has a high resistance when the cap slides behind the patient's head. It is too small !! It is necessary to use a bigger size.

Disposal: After use, dispose the device according to the laws in force in the country of use.

Spes Medica S.r.l.
Via Europa Zona Industriale
84091 Battipaglia (SA) Italy
Tel: +39 0828 614191
Fax: +39 0828 341788
web: www.spesmedica.com
e-mail: info@spesmedica.com

Spes Medica S.r.l.
Via Europa Zona Industriale
84091 Battipaglia (SA) Italy
Tel: +39 0828 614191
Fax: +39 0828 341788
web: www.spesmedica.com
e-mail: info@spesmedica.com

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